






## TYPES OF WAVEFORMS

<u>Type</u>	<u>Clinical Uses</u>	<u>Waveform</u>
<b>MONOPHASIC WAVEFORMS</b>		
Direct Current (DC) or Galvanic	Iontophoresis, Wound Healing, Denervated Tissue	
Pulsed Galvanic or Interrupted DC	Edema Reduction, Wound Healing, Innervated Muscle Contraction	
		
<b>BIPHASIC WAVEFORMS</b>		
Symmetrical AC	Pain Suppression, Innervated Muscle Contraction	
Asymmetrical AC	Pain Suppression, Innervated Muscle Contraction	
<b>TRIPHASIC WAVEFORMS</b>		
Unbalanced Triphasic	Edema Reduction, Pain Suppression	